

Dental Instruments

AI

Impression

Impresión

Eindruck

Quality makes Dentistry Better



Features

No Welding
Easy to Clean

Easy Sterilization/Autoclavable

Sterilizing upto 134C

Adjustable Front Wall

High Quality AISI 304 Stainless Steel

Micro Sand Blast Finish

Also available in Mirror Finished



INOX
STAINLESS STEEL

Impression Trays

New Innovation Impression Trays

ACUMED New and Innovative Impression Trays is the result of ACUMED research laboratories study project aimed of achieving an innovative and economic tool for dental impression taking. Medcave impression trays were designed for non edentulous patients. 5 sizes are available for the upper and lower jaws i.e. (extra small, small, medium, large, extra large).

Acumed Perforated and Non Perforated Impression trays are compatible with all impression materials and available with micro sand blasted stain finish for even more retention.

Acumed addressess all concers raised by plastic and aluminium tray users such as:

1. Not Sufficient handle rigidity.
2. Points of soldering and undercuts hiding place of biological residues.
3. Difficult cleaning processess and ineffective disintection and sterilization procedures.
4. Storing of too many sizes.
5. Need for adhesive products because of not sufficient retension.
6. Patient discomfort.

Disinfecting

Resistant to disinfecting solutions with high quality AISI 304 high grade stainless steel,

Sterilization

Acumed New Impression Trays are fully autoclavable at 134C with no alteration.



Upper

Extra Small	XS	AI-01
Small	S	AI-02
Medium	M	AI-03
Large	L	AI-04
Extra Large	XL	AI-05



Lower

Extra Small	XS	AI-06
Small	S	AI-07
Medium	M	AI-08
Large	L	AI-09
Extra Large	XL	AI-10



Upper

Extra Small	XS	AI-11
Small	S	AI-12
Medium	M	AI-13
Large	L	AI-14
Extra Large	XL	AI-15



Lower

Extra Small	XS	AI-16
Small	S	AI-17
Medium	M	AI-18
Large	L	AI-19
Extra Large	XL	AI-20

Edentulous Impression Trays

Upper - Perforated

Small	S	AI-21
Medium	M	AI-22
Large	L	AI-23

Lower - Perforated

Small	S	AI-24
Medium	M	AI-25
Large	L	AI-26

Upper - Solid

Small	S	AI-27
Medium	M	AI-28
Large	L	AI-29

Lower - Solid

Small	S	AI-30
Medium	M	AI-31
Large	L	AI-32



Impression

Silicon Buttons for Impression Trays


Fit any Impression Tray, easy identification of size.
Assorted pack of 24 buttons, 6 colors, 4 of each.

AI-PACK
(Pack of 24)



ACUMED Metal Impression Trays, the most demanding Trays and first choice of most Dental Schools, Colleges, Universities and Dentists Worldwide. These ACUMED Trays are available in many different sizes, shapes and types including full or partial, Solid or Perforated and Pediatric. When poured and mounted on the articulator, Check Bite Impression produce a pair of contemporary, separable models in true occlusal relationship. High Quality Surgical Grade Stainless Steel is used to manufacture the trays.

Impression Trays Regular, Perforated with Rim Color coded for easy identification of size.



Perforated



Upper



Retention Rim



Lower

Extra Small	XS	AI-33
Small	S	AI-34
Medium	M	AI-35
Large	L	AI-36
Extra Large	XL	AI-37

Extra Small	XS	AI-38
Small	S	AI-39
Medium	M	AI-40
Large	L	AI-41
Extra Large	XL	AI-42



U1 U2 U3 U4 U5



L1 L2 L3 L4 L5

For Ordering Color coding Impression Trays, Please specify.

Impression Trays Regular, Solid with Rim Color coded for easy identification of size.



Solid



Upper



Retention Rim



Lower

Extra Small	XS	AI-43
Small	S	AI-44
Medium	M	AI-45
Large	L	AI-46
Extra Large	XL	AI-47

Extra Small	XS	AI-48
Small	S	AI-49
Medium	M	AI-50
Large	L	AI-51
Extra Large	XL	AI-52



U1 U2 U3 U4 U5



L1 L2 L3 L4 L5

For Ordering Color coding Impression Trays, Please specify.

Stainless Steel Impression Trays for Edentulous



Upper

Extra Small	XS	AI-53
Small	S	AI-54
Medium	M	AI-55
Large	L	AI-56
Extra Large	XL	AI-57



Lower

Extra Small	XS	AI-58
Small	S	AI-59
Medium	M	AI-60
Large	L	AI-61
Extra Large	XL	AI-62



Upper

Extra Small	XS	AI-63
Small	S	AI-64
Medium	M	AI-65
Large	L	AI-66
Extra Large	XL	AI-67



Lower

Extra Small	XS	AI-68
Small	S	AI-69
Medium	M	AI-70
Large	L	AI-71
Extra Large	XL	AI-72

Impression Trays for Partial Impressions

For Children

Upper

Lower



S AI-79	AI-82	S AI-85	AI-88
M AI-80	AI-83	M AI-86	AI-89
L AI-81	AI-84	L AI-87	AI-90



S AI-91	AI-94	S AI-97	AI-100
M AI-92	AI-95	M AI-98	AI-101
L AI-93	AI-96	L AI-99	AI-102

Lateral



Impression Trays

Bite Plane
AI-105



Bite Gouge
AI-106



IMPRESSION TRAYS BITE
CHECA TRAY
AI-107

Classic

Carvers for Silicon Impression Materials

AI-108- 45° Offset



ACUMED Cord Packer # 113

AI-109



ACUMED Cord Packer # 122

AI-110



AI-111



AI-112



ACUMED Cord Packers are also available with finely serrated tips

AI-113- Serrated



AI-114- Smooth





Acumed Logo and Name are Registered and Trade Mark of Acumed



Scan our mobile website to your smartphone for quick reference

Acumed International Co.
Corporate Head
Jinnah Town, Capital Road,
Sialkot - 51310, Pakistan.

Email: sales@acumedco.com

Website: www.acumedco.com